



# CellCor™ CD MSC

## Human Mesenchymal Stem Cells



### Ideal Culture Medium

It is important to select an appropriate medium since it can affect cell properties. CellCor™ CD MSC is the most ideal product to control culture conditions as it consists of chemically defined components.



### Stable Maintenance

CellCor™ CD MSC shows superior proliferation until late passages and stable stemness (Tri-lineage differentiation and expression of MSC markers) in addition to low senescence and genetic stability.



### Universal Use

CellCor™ CD MSC can be used for various tissue-derived MSCs (Adipose tissue, Bone marrow, Umbilical cord, and Wharton's jelly, etc.), and is ideal for exosome research.

## Product Overview

|                                 |                                                                          |
|---------------------------------|--------------------------------------------------------------------------|
| <b>Product</b><br>(Formulation) | <b>CellCor™ CD MSC</b><br>Chemically Defined (Recombinants / Synthetics) |
| <b>Catalog</b>                  | YSP002                                                                   |
| <b>Application</b>              | Cell culture and expansion                                               |
| <b>Cell Type</b>                | Human Mesenchymal Stem Cells (MSCs)                                      |
| <b>Storage</b>                  | Under -20°C (Expiry date on label)                                       |
| <b>Supplement</b>               | No supplement and coating reagent                                        |
| <b>Size</b>                     | 500 mL (Custom packaging available)                                      |
| <b>Area</b>                     | For Research and Further Manufacturing Use                               |



# Why CellCor™ CD MSC?



## Defined & Regulatory

Regulatory-friendly as it does not contain animal or human-derived components such as sera or hPL.

## Stemness

Maintains stable stemness until late passages.  
Genetic stability at the nuclear and chromosomal level in cell culture.

## Consistency & Quality

Chemically defined products provide consistent products even in mass production.

## Yield & Productivity

Excellent cell productivity with superior proliferation.


## Reproducibility

Able to obtain highly reproducible experimental results.


# MSC Stemness

Superior & Standard Stem Cell Characterization


## Cell Morphology



## Proliferation



## Colony Formation




### Figure 1. Superior Proliferation

AdMSCs cultured with CellCor™ CD MSC shows typical cell morphology, superior proliferation, and superior colony formation. (CellCor™ CD MSC : Chemically Defined Medium, FBS : FBS Containing Medium)

\* Reference : Humana Press, 2008. 83-91


## Tri-differentiation

(Passage 7)



## MSC Marker

(Passage 7)




### Figure 2. Stable MSC characterization

AdMSCs cultured with CellCor™ CD MSC maintains stemness : tri-lineage differentiation (adipogenesis, chondrogenesis, and osteogenesis) and MSC specific marker. (CD14, CD34, CD45, HLA-DR, CD73, CD90, and CD105) (Passage 7)

\* Reference : Cytotherapy, 21 (10), 1019-1024


## Senescence (Staining)

(Passage 5)



## Senescence Marker (ODZ2)

(Passage 5)




### Figure 3. Low Senescence

AdMSCs cultured with CellCor™ CD MSC shows less b-galactosidase staining indicative of cellular senescence. (Passage 5) (CellCor™ CD MSC : Chemically Defined Medium, FBS: FBS Containing Medium)

\* Reference : Cells. 2021 Jun; 10 (6) : 1301.

# Genetic Stability


Minimize loss of genetic information



### Figure 4. CBMN assay

AdMSCs cultured with CellCor™ CD MSC shows genetic stability at the nuclear level in cell culture. (Passage 6)  
(CellCor™ CD MSC : Chemically Defined Medium, FBS: FBS Containing Medium)

\* Reference : Cytotherapy, 21 (10), 1019–1024.




### Figure 5. Karyotype analysis

AdMSCs cultured with CellCor™ CD MSC shows genetic stability at the chromosomal level in cell culture. (Passage 5)  
(CellCor™ CD MSC : Chemically Defined Medium, FBS: FBS Containing Medium)

\* Reference : Cytotherapy 18.3 (2016): 336–343.

# Defined & Consistency

## Unique features of Chemically Defined medium




**Figure 1. Superior Proliferation**

CellCor™ CD MSC is a chemically defined medium that shows it does not contain serum or animal-derived components. (Passage 5, Neu5Gc: Animal-derived marker, BDKRB1: Human-derived marker). (CellCor™ CD MSC: Chemically Defined Medium, FBS: FBS Containing Medium, Human Serum: Human Serum Containing Medium)

\* Reference : Tissue Eng Part A, 2010, 16(4): p. 1143-55 / international journal of medical sciences 16.8 (2019): 1102.

## Lot to Lot Variation

(Lot# N = 5)




**Figure 2. Consistency**

CellCor™ CD MSC, Chemically defined medium, shows little variation as it does not contain animal- or human-derived components. It can improve the reproducibility of the experiment. (compared to serum-free medium; CellCor™ SFD MSC)

# General & Universal

Applicable to various cell types




**Figure 1.** Characteristics of MSC from various tissue source

CellCor™ CD MSC could be cultured in various tissue- and cell-derived sources. (Adipose tissue, Umbilical cord, Bone marrow, Wharton's Jelly, iPSC) and shows to maintain stemness. (CellCor™ CD MSC : Chemically Defined Medium, FBS : FBS Containing Medium, Human Serum : Human Serum Containing Medium)





# Chemically Defined vs Serum-free Media

## Cell Proliferation




## Senescence Staining

(Passage 5)




## MSC Marker

(Passage 5)




## DEG

(Differentially Expressed Genes) (Passage 5)




## Tri-differentiation

(Passage 5)



## Genetic Stability

(Passage 5)



**Figure 1.** Compare MSC characterization to Serum-free Medium

Comparing with serum-free medium (Commercial product) AdMSCs cultured with CellCor CD MSC show superior cell proliferation, stable stemness, low senescence. Also, it shows low expression of immune- and inflammatory-related genes, and cell culture of genetic stability. (Passage 5, Ag-1 : IGF1, Ag-2 : Endoglin, Ag-3 : SOD2, If-1 : TNFRSF11B, If-2 : CSF1, Im-1 : CXCL12, Im-2 : PTX3)



# CellCor™ CD MSC Applications





CellCor™ is a chemically defined medium specifically for cell culture.  
With CellCor™ CD MSC as a start, we are constantly developing ways to expand  
our product line to various cell specific chemically defined medium.

## Technical Support

---

E. [ts@xcell.co.kr](mailto:ts@xcell.co.kr)

## Website Resources

---

Material Safety Data Sheets (MSDS)  
Certificate of Analysis (CoA)  
Product Information Sheets (PIS)  
Product Manual Video  
FAQs

## Xcell Therapeutics

---

### Head Office (SEOUL)

06188  
Dongwon Bldg. 6F, 333, Yeongdong-daero,  
Gangnam-gu, Seoul, Korea

T. +82-2-863-1331

F. +82-2-863-0832

E. [info@xcell.co.kr](mailto:info@xcell.co.kr)

### GMP Facility (YONGIN)

16954  
#1908, 13, Heungdeok 1-ro, Giheung-gu,  
Yongin-si, Gyeonggi-do, Korea

T. +82-31-214-1088

F. +82-31-660-7058

E. [gmp@xcell.co.kr](mailto:gmp@xcell.co.kr)

