

Compatibility Chart 2026-06

GenFollower Biotech CO., LTD.

ADD: Eastern industrial park, Baiguan street, Shangyu district,
Shaoxing city, Zhejiang province, China
Web: www.GenFollower.com



FILTER PIPETTE TIPS

CAT.#	Eppendorf research	Eppendorf research plus	Eppendorf reference	Gilson pipetman	Finpipette F3	Finpipette	Biohit	CAPP	SOCOREX	Brand Tripette
E-FTB10S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-FTB10S-L	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-FTR10S-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-FTR10S-L-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-FTB10GL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-FTB10GL-L	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-FTR10GL-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-FTR10GL-L-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-FTB20	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTB20-L	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTR20-S	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTR20-L-S	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTB50	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTB50-L	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTR50-S	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTR50-L-S	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTB100	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	10µL - 100µL, 20µL - 200µL	10µL - 100µL, 20µL - 200µL	P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTB100-L	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	10µL - 100µL, 20µL - 200µL	10µL - 100µL, 20µL - 200µL	P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL

Remark: Considering the size of different brands of pipettes might be changed year to year, so please kindly test your pipettes with GenFollower tips before purchasing large quantity.

Compatibility Chart 2026-06

GenFollower Biotech CO., LTD.
ADD: Eastern industrial park, Baiguan street, Shangyu district,
Shaoxing city, Zhejiang province, China
Web: www.GenFollower.com



FILTER PIPETTE TIPS

CAT.#	Eppendorf research	Eppendorf research plus	Eppendorf reference	Gilson pipetman	Finpipette F3	Finpipette	Biohit	CAPP	SOCOREX	Brand Tripette
E-FTR100-S	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	10µL - 100µL, 20µL - 200µL	10µL - 100µL, 20µL - 200µL	P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTR100-L-S	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	10µL - 100µL, 20µL - 200µL	10µL - 100µL, 20µL - 200µL	P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTB200	30µL - 300µL	30µL - 300µL	30µL - 300µL	P200	20µL-200µL, 50µL-300µL	20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTB200-L	30µL - 300µL	30µL - 300µL	30µL - 300µL	P200	20µL-200µL, 50µL-300µL	20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTR200-S	30µL - 300µL	30µL - 300µL	30µL - 300µL	P200	20µL-200µL, 50µL-300µL	20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTR200-L-S	30µL - 300µL	30µL - 300µL	30µL - 300µL	P200	20µL-200µL, 50µL-300µL	20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTB300	30µL - 300µL	30µL - 300µL	30µL - 300µL	P200	20µL-200µL, 50µL-300µL	20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTB300-L	30µL - 300µL	30µL - 300µL	30µL - 300µL	P200	20µL-200µL, 50µL-300µL	20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTR300-S	30µL - 300µL	30µL - 300µL	30µL - 300µL	P200	20µL-200µL, 50µL-300µL	20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTR300-L-S	30µL - 300µL	30µL - 300µL	30µL - 300µL	P200	20µL-200µL, 50µL-300µL	20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-FTB1000	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-FTB1000-L	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-FTR1000-S	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-FTR1000-L-S	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-FTB1250	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-FTB1250-L	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-FTR1250-S	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-FTR1250-L-S	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL

Remark: Considering the size of different brands of pipettes might be changed year to year, so please kindly test your pipettes with GenFollower tips before purchasing large quantity.

Compatibility Chart 2026-06

GenFollower Biotech CO., LTD.
ADD: Eastern industrial park, Baiguan street, Shangyu district,
Shaoxing city, Zhejiang province, China
Web: www.GenFollower.com



FILTER PIPETTE TIPS

CAT.#	Eppendorf research	Eppendorf research plus	Eppendorf reference	Gilson pipetman	Finpipette F3	Finpipette	Biohit	CAPP	SOCOREX	Brand Tripette
E-DFTR10GL-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-DFTR100-S	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	10µL - 100µL, 20µL - 200µL	10µL - 100µL, 20µL - 200µL	P100, P200	5µL-50µL, 20µL-200µL, 50µL-300µL	5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-DFTR200-S	30µL - 300µL	30µL - 300µL	30µL - 300µL	P200	20µL-200µL, 50µL-300µL	20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-DFTR300-S	30µL - 300µL	30µL - 300µL	30µL - 300µL	P200	20µL-200µL, 50µL-300µL	20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL

Remark: Considering the size of different brands of pipettes might be changed year to year, so please kindly test your pipettes with GenFollower tips before purchasing large quantity.

Compatibility Chart 2026-06

GenFollower Biotech CO., LTD.
ADD: Eastern industrial park, Baiguan street, Shangyu district,
Shaoxing city, Zhejiang province, China
Web: www.GenFollower.com



PIPETTE TIPS WITHOUT FILTER

CAT.#	Eppendorf research	Eppendorf research plus	Eppendorf reference	Gilson pipetman	Finpipette F3	Finpipette	Biohit	CAPP	SOCOREX	Brand Tripette
E-TB10S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-TB10S-L	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-TR10S-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-TR10S-L-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-TB10GL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-TB10GL-L	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-TR10GL-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-TR10GL-L-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-RS10GL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-RS10GL-L	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-RS10GL-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-RS10GL-L-S	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	0.1µL - 2.5µL, 0.5µL - 10µL	P10	/	/	Proline Plus 0.5-10uL	0.5-10uL	1-10µL	0.5-10uL
E-TB200	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	0.5µL-10µL, 5µL-50µL, 20µL-200µL, 50µL-300µL	0.5µL-10µL, 5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-TB200-L	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	0.5µL-10µL, 5µL-50µL, 20µL-200µL, 50µL-300µL	0.5µL-10µL, 5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-TR200-S	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	0.5µL-10µL, 5µL-50µL, 20µL-200µL, 50µL-300µL	0.5µL-10µL, 5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL
E-TR200-L-S	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	2µL - 20µL, 10µL - 100µL, 20µL - 200µL	P20, P100, P200	0.5µL-10µL, 5µL-50µL, 20µL-200µL, 50µL-300µL	0.5µL-10µL, 5µL-50µL, 20µL-200µL, 50µL-300µL	M100/M200	5uL-50uL, 20uL-200uL	50µL-200µL	20µL-200µL, 30µL-300µL

Remark: Considering the size of different brands of pipettes might be changed year to year, so please kindly test your pipettes with GenFollower tips before purchasing large quantity.

Compatibility Chart 2026-06

GenFollower Biotech CO., LTD.

ADD: Eastern industrial park, Baiguan street, Shangyu district,
Shaoxing city, Zhejiang province, China
Web: www.GenFollower.com



PIPETTE TIPS WITHOUT FILTER

CAT.#	Eppendorf research	Eppendorf research plus	Eppendorf reference	Gilson pipetman	Finpipette F3	Finpipette	Biohit	CAPP	SOCOREX	Brand Tripette
E-RS200	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-RS200-L	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-RS200-S	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-RS200-L-S	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-TB300	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-TB300-L	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-TR300-S	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-TR300-L-S	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-RS300	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-RS300-L	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-RS300-S	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-RS300-L-S	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	2μL - 20μL, 10μL - 100μL, 20μL - 200μL	P20, P100, P200	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	0.5μL-10μL, 5μL-50μL, 20μL-200μL, 50μL-300μL	M100/M200	5uL-50uL, 20uL-200uL	50μL-200μL	20μL-200μL, 30μL-300μL
E-TB1000	100μL - 1000μL	100μL - 1000μL	100μL - 1000μL	P1000	100-1000μL/200-1000μL	100-1000μL/200-1000μL	/	100-1000μL	100-1000μL	100-1000μL
E-TB1000-L	100μL - 1000μL	100μL - 1000μL	100μL - 1000μL	P1000	100-1000μL/200-1000μL	100-1000μL/200-1000μL	/	100-1000μL	100-1000μL	100-1000μL
E-TR1000-S	100μL - 1000μL	100μL - 1000μL	100μL - 1000μL	P1000	100-1000μL/200-1000μL	100-1000μL/200-1000μL	/	100-1000μL	100-1000μL	100-1000μL
E-TR1000-L-S	100μL - 1000μL	100μL - 1000μL	100μL - 1000μL	P1000	100-1000μL/200-1000μL	100-1000μL/200-1000μL	/	100-1000μL	100-1000μL	100-1000μL

Remark: Considering the size of different brands of pipettes might be changed year to year, so please kindly test your pipettes with GenFollower tips before purchasing large quantity.

Compatibility Chart 2026-06

GenFollower Biotech CO., LTD.
ADD: Eastern industrial park, Baiguan street, Shangyu district,
Shaoxing city, Zhejiang province, China
Web: www.GenFollower.com



PIPETTE TIPS WITHOUT FILTER

CAT.#	Eppendorf research	Eppendorf research plus	Eppendorf reference	Gilson pipetman	Finpipette F3	Finpipette	Biohit	CAPP	SOCOREX	Brand Tripette
E-RS1000	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-RS1000-L	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-RS1000-S	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-RS1000-L-S	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-TB1250	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-TB1250-L	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-TR1250-S	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-TR1250-L-S	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-RS1250	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-RS1250-L	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-RS1250-S	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL
E-RS1250-L-S	100µL - 1000µL	100µL - 1000µL	100µL - 1000µL	P1000	100-1000µL/200-1000µL	100-1000µL/200-1000µL	/	100-1000µL	100-1000µL	100-1000µL

Remark: Considering the size of different brands of pipettes might be changed year to year, so please kindly test your pipettes with GenFollower tips before purchasing large quantity.