

## Weighing Dish

(#99010)

### 1. Description

SPL weighing dishes with channels for material flow provide better reagent loss through antistatic performance.

### 2. Features

- ❖ Manufactured from modified polystyrene
- ❖ Material flow channel towards the edge
- ❖ Antistatic

### 3. Product Range

Cat. No.	Material	Color	Internal Dimension w x l(mm)	Volume(ml)	Sterile	Packaging
99010	PS	White	77 x 77	100.00	-	125 / 250

### 5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

### 6. Product Photographs




7. Specifications & Technical Drawing

- Specifications

	Dimensions(External)	Material (Resin)	Packaging
	99010		
Length:	89 mm	Antistatic treated polystyrene (PS)	125 / 250
Width:	89 mm		
Height	25 mm		
Working volume:	100 ml		

- Technical Drawing

Cat. No. 99010



[www.spllifesciences.com](http://www.spllifesciences.com)

For technical assistance, contact SPL R&D Center at:  
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [spl@ispl.co.kr](mailto:spl@ispl.co.kr)  
 To place an order, contact your local distributor or  
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [business@ispl.co.kr](mailto:business@ispl.co.kr)

SPL-TDS-RES

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.  
 © SPL Life Sciences Co., Ltd., Version 2, July 20<sup>th</sup>, 2016

