

Certification of Analysis

Product Name: CD20-PE

Clone: HI47

Batch Number: 202209070

Quantity: 100T

Production date: 20221013

Expiration Date: 20240413

Buffer and Preservative: 0.01M PBS (pH 7.2) with 0.1% proclin300 and 0.2% BSA

Method of Test: Direct immunofluorescence staining for flow cytometric analysis

Suggested working titer: 20 μ l/1 \times 10⁶ cells or 100 μ l peripheral blood

Results: This batch of products reacts with normal human peripheral blood lymphocytes (18.9%), is negative in peripheral blood monocytes and granulocytes.

Conclusion: This batch of products conforms to the specificity of anti-human CD20 antibody.

Technician: Meng yuan Han Date: 10/13/2022

Authorized for release: Ying Mo Date: 10/13/2022

www.quantobio.com

E-mail: info@quantobio.com

Tel: +86 10 59773 899

Fax: +86 10 59773 998

For research use only. Not for diagnostic use. Not for resale. Quantobio will not be held responsible for patent infringement or other violations that may occur with the use of our products.