

Certification of Analysis

Product Name: CD34-PE

Clone: 4H11

Batch Number: 202212031

Quantity: 100T

Production date: 20230103

Expiration Date: 20240703

Buffer and Preservative: 0.01M PBS (pH 7.2) with 0.1% proclin300 and 0.2% BSA

Method of Test: Direct immunofluorescence staining for flow cytometric analysis

Suggested working titer: 20 μ l/1 \times 10⁶ cells or 100 μ l peripheral blood

Results: This batch of products reacts with KG1a cell lines (100%), is negative in normal human peripheral blood lymphocytes, granulocytes and monocytes.

Conclusion: This batch of products conforms to the specificity of anti-human CD34 antibody.

Technician: Mengyuan Han Date: 01/03/2023

Authorized for release: Ying Ma Date: 01/03/2023

www.quantobio.com

E-mail: info@quantobio.com

Tel: +86 10 59773 899

Fax: +86 10 59773 998

For research use only. Not for diagnostic use. Not for resale. Quantobio will not be held responsible for patent infringement or other violations that may occur with the use of our products.