

Certification of Analysis

Product Name: CD45-FITC

Clone: HI30

Batch Number: 202209077

Quantity: 100T

Production date: 20221017

Expiration Date: 20240417

Buffer and Preservative: 0.01M PBS (pH 7.2) with 0.1% proclin300 and 0.2% BSA

Method of Test: Direct immunofluorescence staining for flow cytometric analysis

Suggested working titer: 20 μ l/1 \times 10⁶ cells or 100 μ l peripheral blood

Results: This batch of products reacts with normal human peripheral blood lymphocytes (100%), granulocytes (100%), monocytes (100%).

Conclusion: This batch of products conforms to the specificity of anti-human CD45 antibody.

Technician: Mengyuan Han Date: 10/17/2022

Authorized for release: Ying Chen Date: 10/17/2022

www.quantobio.com

E-mail: info@quantobio.com

Tel: +86 10 59773 899

Fax: +86 10 59773 998

For research use only. Not for diagnostic use. Not for resale. Quantobio will not be held responsible for patent infringement or other violations that may occur with the use of our products.