

Certification of Analysis

Product Name: Purified Mouse Anti-human CD9 mAb

Clone: HI9a

Batch Number: P202106019

Amount: 0.1mg

Concentration: 1.0mg/ml

Production date: 20210616

Expiration Date: 20240616

Purity: $\geq 90\%$ by SDS-PAGE

Buffer and Preservative: 0.01M PBS pH7.2 with 0.1% NaN₃

Method of Purification: Affinity chromatography

Method of Test: Indirect immunofluorescence staining for flow cytometric analysis

Suggested working titer: 0.25 μ g/ 1×10^6 cells or 100 μ l peripheral blood

Results: This batch of products reacts with normal human platelets (75.14%).

Conclusion: This batch of products conforms to the specificity of anti-human CD9 antibody.

Technician: Dishan Date: 06/16/2021

Authorized for release: Ying Ma Date: 06/16/2021

www.quantobio.com

E-mail: info@quantobio.com

Tel: +86 10 59773 899

Fax: +86 10 59773 998

For research use only. Not for diagnostic use. Not for resale. Quantobio will not be held responsible for patent infringement or other violations that may occur with the use of our products.